

cowie® PTFE Stirring Magnets


HOME

LABWARE
INDEX

For profile drawing see:
Plain
Pivot Ring
Triangle


Giant

For use in very large containers



Oval

For use in round bottom flasks



Tapered

An effective stirrer with similar action to Oval/Piv of Ring type

Reference	L (mm) x D (mm)	Reference	L (mm) x D (mm)	Reference	L (mm) x D (mm)
Plain					
001.257	57 x 27	001.610	10 x 5	001.1910	10 x 4
001.2108	108 x 27	001.615	15 x 6	001.1915	15 x 5
001.2159	159 x 27	001.620	20 x 10	001.1920	20 x 7
		001.625.10	25 x 10	001.1925	25 x 8
		001.625	25 x 12	001.1930	30 x 8
		001.630.10	30 x 10	001.1935	35 x 8
		001.630	30 x 16	001.1940	40 x 8
		001.635.13	35 x 13	001.1945	45 x 8
		001.635	35 x 16	001.1950	50 x 8
		001.640.13	40 x 13	001.1955	55 x 8
		001.640	40 x 20	001.1960	60 x 8
		001.650.17	50 x 17	001.1970	70 x 10
		001.650	50 x 20	001.1980	80 x 10
		001.664	64 x 20		
		001.670	70 x 20		
		001.670.25	70 x 25		
		001.670.27	70 x 27		
		001.6100	100 x 30		
		001.6150	150 x 40		
Pivot Ring					
001.357	57 x 24				
001.3108	108 x 24				
001.3159	159 x 24				
Triangle					
001.480	80 x 18				
001.4110	108 x 36				
001.4136	136 x 36				

NB: The dimensions of stirrer bars are nominal and approximate to the stated value by $\pm 5\%$ for length and $\pm 10\%$ for diameter

cowie®